

Crimpcaps with septa (11 mm)



					
Article number	150002	150004	150005	150006	150007
Description	Crimpcap 11mm, blank 1.0 mm septa, red NR / butylrubber / PTFE	Crimpcap 11mm, red 1.0 mm septa, red NR / butylrubber / PTFE	Crimpcap 11mm, green 1.0 mm septa, red NR / butylrubber / PTFE	Crimpcap 11mm, blue 1.0 mm septa, red NR / butylrubber / PTFE	Crimpcap 11mm, gold 1.0 mm septa, red NR / butylrubber / PTFE
Packaging	1000 pcs / bag	1000 pcs / bag	1000 pcs / bag	1000 pcs / bag	1000 pcs / bag

					
Article number	151130	151133	151131	151132	151134
Description	Crimpcap 11mm, blank 1.2 mm septa, white / red Sil / PTFE	Crimpcap 11mm, red 1.2 mm septa, white / red Sil / PTFE	Crimpcap 11mm, green 1.2 mm septa, white / red Sil / PTFE	Crimpcap 11mm, blue 1.2 mm septa, white / red Sil / PTFE	Crimpcap 11mm, gold 1.2 mm septa, white / red Sil / PTFE
Packaging	1000 pcs / bag	1000 pcs / bag	1000 pcs / bag	1000 pcs / bag	1000 pcs / bag

Crimpcaps with septa (11 mm)

					
Article number	151135	151140	151166		
Description	Crimpcap 11 mm, black 1.2 mm septa, white / red Sil / PTFE	Crimpcap 11 mm, blank with rill-rand PTFE liner, white, 0.25 mm	Crimpcap 11 mm, blank with rill-rand PTFE liner, white, 0.25 mm, with TPF O-ring		
Packaging	1000 pcs / bag	1000 pcs / bag	1000 pcs / bag		