









Crimptop vials (11 mm)



				
Article number	110006	110018	110019	110007
Description	1.5 ml crimptop vial 32 * 11.6 mm glass, clear flat bottom	0.3 ml crimptop vial 32 * 11.6 mm polypropylene, clear flat bottom	1.5 ml crimptop vial 32 * 11.6 mm glass, clear wide mouth flat bottom	1.5 ml crimptop vial 32 * 11.6 mm glass, clear wide mouth flat bottom white 'write on' patch
Packaging	100 pcs / box	100 pcs / box	100 pcs / box	100 pcs / box
				
Article number	110048	110004	110022	110046
Description	0.2 ml crimptop vial 32 * 11.6 mm glass, clear 'write on' patch integrated insert base bonded	1.5 ml crimptop vial 32 * 11.6 mm glass, amber narrow mouth (5 mm) flat bottom	0.2 ml crimptop vial 32 * 12 mm glass, amber flat bottom insert integrated top bonded	1.5 ml crimptop vial 32 * 11.6 mm glass, amber wide mouth flat bottom white 'write on' patch
Packaging	100 pcs / box	100 pcs / box	100 pcs / box	100 pcs / box

Crimptop vials (11 mm)

				
Article number	110047	110042	110043	110044
Description	0.2 ml crimptop vial 32 * 11.6 mm glass, amber 'write on' patch integrated insert, base bonded	1.5 ml crimptop vial 32 * 11.6 mm glass, clear wide mouth flat bottom red 'write on' patch	1.5 ml crimptop vial 32 * 11.6 mm glass, clear wide mouth flat bottom blue 'write on' patch	1.5 ml crimptop vial 32 * 11.6 mm glass, clear wide mouth flat bottom yellow 'write on' patch
Packaging	100 pcs / box	100 pcs / box	100 pcs / box	100 pcs / box
				
Article number	110045	110021		
Description	1.5 ml crimptop vial 32 * 11.6 mm glass, clear wide mouth flat bottom green 'write on' patch	1.1 ml crimptop vial 32 * 11.6 mm glass, clear narrow mouth conical		
Packaging	100 pcs / box	100 pcs / box		