





Screwcap vials (13 mm)



				
Article number	120031	120032	120033	120034
Description	4 ml screwcap vial 13-425 45*15 mm glass, clear flat bottom	4 ml screwcap vial 13-425 45*15 mm glass, amber flat bottom	4 ml screwcap vial 13-425 45*15 mm glass, clear flat bottom write on patch	4 ml screwcap vial 13-425 45*15 mm glass, amber flat bottom write on patch
Packaging	100 pcs / package	100 pcs / package	100 pcs / package	100 pcs / package