


# 13 mm screwcaps



				
<b>Article number</b>	970287	151120		
<b>Description</b>	Screwcap for 4 ml vial black, polypropylene Sil/PTFE septa transp./white	Screwcap for 4 ml vial black, polypropylene Sil/PTFE septa white/blue pre-slit		
<b>Packaging</b>	100 pcs / bag	1000 pcs / bag		
				
<b>Article number</b>	151045	900486	900988	151046
<b>Description</b>	Screwcap for 4 ml vial black, polypropylene with open top	Screwcap for vial 13-425 black, butyl/PTFE closed	Screwcap for vial 13-425 black, PE/PTFE closed	Screwcap for 4 ml vial black, bakelite closed
<b>Packaging</b>	1000 pcs / bag	1000 pcs / bag	1000 pcs / bag	1000 pcs / bag