

Vial & cap 13 mm

screwcap



120031



120032



120033



120034

Article	Description	Packaging
120031	Vial, screwcap, 13-425, 4 ml, 45*15 mm, glass, clear, flat bottom	100
120032	Vial, screwcap, 13-425, 4 ml, 45*15 mm, glass, amber, flat bottom	100
120033	Vial, screwcap, 13-425, 4 ml, 45*15 mm, glass, clear, flat bottom, writing patch	100
120034	Vial, screwcap, 13-425, 4 ml, 45*15 mm, glass, amber, flat bottom, writing patch	100



151045



970287



151046



900486

Article	Description	Packaging
151045	Cap, screw, for vial, 4 ml, 13-425, black, polypropylene, with open top	100
970287	Cap, screw, for vial, 4 ml, 13-425, black, polypropylene, Sil/PTFE	100
151046	Cap, screw, for vial, 4 ml, 13-425, black, bakelite, closed	100
900486	Cap, screw, for vial, 13-425, black, bakelite, closed, butyl/PTFE	100