




# Screwcap vials and caps GL18 (18 mm)



				
<b>Product number</b>	120037	120046	120059	970266
<b>Description</b>	9 ml screwcap vial, GL18 (18 mm) 40 * 22 mm glass, clear flat bottom	16 ml screwcap vial, GL18 (18 mm) 71 * 20.6 mm glass, clear	20 ml screwcap vial, GL18 (18 mm) 75.5 * 22.5 mm glass, clear	16 ml screwcap vial, GL18 (18 mm), assembled 71 * 20.6 mm glass, clear screwcap foam/PTFE septa
<b>Packaging</b>	186 pieces/package	100 pieces/package	100 pieces/package	100 pieces/package
				
<b>Product number</b>	151100	151142	151143	151101
<b>Description</b>	screwcap for 9 ml GL18 (18 mm) vial white, polypropylene foam septa	screwcap for GL18 (18 mm) vial black, polypropylene with injectionhole	septa 16.5 * 2 mm PTFE/EPE brown/white	septa 17 * 0.2 mm PTFE white
<b>Packaging</b>	individually	7000 pieces/package	1000 pieces/package	1000 pieces/package