











Crimptop vials (20 mm)



				
Product number	112001	112010	112020	112026
Description	5 ml headspace vial 38 * 22 mm glass, clear flat bottom flat top	10 ml headspace vial 45 * 22 mm glass, clear flat bottom flat top	20 ml headspace vial 76 * 23 mm glass, clear flat bottom flat top	20 ml headspace vial 75.5 * 22.5 mm glass, clear round bottom flat top
Packaging	100 pcs / package	100 pcs / package	100 pcs / package	100 pcs / package

			
Product number	112023	112029	112040
Description	20 ml headspace vial 75.5 * 23 mm glass, clear round bottom bevelled top	20 ml headspace vial 75.5 * 22.5 mm glass, clear round bottom bevelled top write on patch	22 ml headspace vial 85 * 22.5 mm glass, clear flat bottom top = 10 mm
Packaging	100 pcs / package	100 pcs / package	100 pcs / package

Crimptop vials (20 mm)

			
Product number	112012	112027	112030
Description	10 ml headspace vial 45 * 22 mm glass, amber flat bottom flat top	20 ml headspace vial 75.5 * 22.5 mm glass, amber round bottom flat top	20 ml headspace vial 75.5 * 22.5 mm glass, amber round bottom bevelled top
Packaging	100 pcs / package	100 pcs / package	100 pcs / package