

Vials 20 mm headspace



112001



112004



112010



112012



112011



112014



112015



900241

Article	Description	Packaging
112001	Vial, headspace, 5 ml, 38*22 mm, glass, clear, flat bottom	100
112004	Vial, headspace, 5 ml, 38*20 mm, glass, clear, flat bottom	100
112010	Vial, headspace, 10 ml, 45*22 mm, glass, clear, flat bottom	100
112012	Vial, headspace, 10 ml, 45*22 mm, glass, amber, flat bottom	100
112011	Vial, headspace, 10 ml, 45*22 mm, glass, clear, flat bottom, 6 mm flange	100
112014	Vial, headspace, 10 ml, 46*22,5 mm, glass, clear, rounded bottom, DIN crimpneck	100
112015	Vial, headspace, 10 ml, 46*22 mm, glass, transparent, bevel top, RB	100
900241	Vial, headspace, 10 ml, 45*22, glass, clear, cleaned	100



112026



112027



112023



112030



112029



112020



112040

Article	Description	Packaging
112026	Vial, headspace, 20 ml, 75.5*22.5 mm, glass, clear, round bottom, flat top	100
112027	Vial, headspace, 20 ml, 75,5*22,5 mm, glass, amber, round bottom, flat top	100
112023	Vial, headspace, 20 ml, 75,5*22,5 mm, glass, clear, round bottom, bevel top	100
112030	Vial, headspace, 20 ml, 75,5*22,5 mm, glass, amber, round bottom, bevel top	100
112029	Vial, headspace, 20 ml, 75.5*22.5 mm, glass, clear, writing patch, round bottom, bevel top	100
112020	Vial, headspace, 20 ml, 76*23 mm, glass, clear, flat bottom, flat top	100
112040	Vial, headspace, 22 ml, 85*22.5 mm, glass, clear, flat bottom, top = 10 mm	100