

Vials (EPA) 24 mm

Screwcap



Article	Description	Packaging
120055	Vial, EPA, 24-400, 20 ml, 57*27.5 mm, glass, clear, flat bottom	100
120051	Vial, EPA, 24-400, 20 ml, 57*27.5 mm, glass, amber, flat bottom	100
120091	Vial, EPA, 24-400, 30 ml, 72.5*27.5 mm, glass, clear, flat bottom	100
120092	Vial, EPA, 24-400, 30 ml, 72.5*27.5 mm, glass, amber, flat bottom	100
120061	Vial, EPA, 24-400, 40 ml, 95*27.5 mm, glass, clear, flat bottom	100
120062	Vial, EPA, 24-400, 40 ml, 95*27.5 mm, glass, amber, flat bottom	100
120064	Vial, EPA, 24-400, 40 ml, 95*27.5 mm, glass, clear, flat bottom, graduated at 25 ml	100
120070	Vial, EPA, 24-400, 40 ml, 82*29 mm, glass, clear, flat bottom, for mobile GC	100
120071	Vial, EPA, 24-400, 60 ml, 140*27.5 mm, glass, clear, flat bottom	100
120072	Vial, EPA, 24-400, 60 ml, 140*27.5 mm, glass, amber, flat bottom	100
120073	Vial, EPA, 24-400, 60 ml, 140*27.5 mm, glass, clear, flat bottom, with graduation	100



Article	Description	Assembled	Packaging
970214	Vial, EPA, 24-400, assembled, 40 ml, glass, clear, Sil/PTFE, transp./ white		100
970419	Vial, EPA, 24-400, assembled, 40 ml, glass, amber, Sil/PTFE, beige/wit		100
970289	Vial, EPA, 24-400, assembled, 40 ml, glass, amber, Sil/PTFE, beige/white, green cap		100
970290	Vial, EPA, 24-400, assembled, 40 ml, glass, amber, Sil/PTFE, beige/white, red cap		100
970300	Vial, EPA, 24-400, assembled, 40 ml, glass, clear, foam/PTFE, white/brown		100
970262	Vial, EPA, 24-200, assembled, 60 ml, glass, amber, Sil/PTFE beige/wit		100
970285	Vial, EPA, 24-400, assembled, 40 ml, glass, amber, foam/PTFE, white/brown		100