









# Screwcap vials (8 mm)



				
<b>Article number</b>	120003	120002	120006	120004
<b>Description</b>	1.5 ml screwcap vial 32*11.6 mm glass, clear narrow mouth, 8-425 flat bottom	1.5 ml screwcap vial 32*11.6 mm glass, amber narrow mouth, 8-425 flat bottom	1.5 ml screwcap vial 32*11.6 mm glass, clear narrow mouth, 8-425 flat bottom write on patch	1.5 ml screwcap vial 32*11.6 mm glass, amber narrow mouth, 8-425 flat bottom write on patch
<b>Packaging</b>	100 pcs / box	100 pcs / box	100 pcs / box	100 pcs / box
				
<b>Article number</b>	120021	120022	120011	120012
<b>Description</b>	1.1 ml screwcap vial 32*12 mm glass, clear 8-425 conical	1.1 ml screwcap vial 32*12 mm glass, amber 8-425 conical	0.2 ml screwcap vial 32*12 mm TPX, clear 8-425 flat bottom	0.2 ml screwcap vial 32*12 mm PP, natural 8-425 flat bottom with graduation
<b>Packaging</b>	100 pcs / box	100 pcs / box	1000 pcs / box	1000 pcs / box