

Vial 8 mm screwcap



120003



120002



120006



120017



970186

Article	Description	Packaging
120003	Vial, screwcap, 8-425, 1.5 ml, 32*11.6 mm, glass, clear, flat bottom, NM	100
120002	Vial, screwcap, 8-425, 1.5 ml, 32*11.6 mm, glass, amber, flat bottom, NM	100
120006	Vial, screwcap, 8-425, 1.5 ml, 32*11.6 mm, glass, clear, flat bottom, writing patch, NM	100
120017	Vial, screwcap, 8-425, 1.5 ml, 32*11.6 mm, glass, amber, flat bottom, writing patch	100
970186	Vial, screwcap, 8-425, cap, liner, 2 ml, Sil/PTFE, assembled	100



151156



151157



151158



151159



970247



151018

Article	Description	Packaging
151156	Cap, screw, for vial, 2 ml, 8-425, red, polypropylene, Sil/PTFE, transparent/red	1000
151157	Cap, screw, for vial, 2 ml, 8-425, white, polypropylene, Sil/PTFE, transparent/red	1000
151158	Cap, screw, for vial, 2 ml, 8-425, green, polypropylene, Sil/PTFE, transparent/red	1000
151159	Cap, screw, for vial, 2 ml, 8-425, black, polypropylene, Sil/PTFE transparent/red	1000
970247	Cap, screw, for vial, 2 ml, 8-425, black, polypropylene, Sil/PTFE, red/white, pre-slit 4 mm	1000
151018	Cap, screw, for vial, 2 ml, 8-425, black, polypropylene, closed cap	1000