








Screwcap vials (9 mm)



				
Article number	120009	120016	120014	120015
Description	1.5 ml screwcap vial 32*11.6 mm glass, clear narrow mouth flat bottom	1.5 ml screwcap vial 32*11.6 mm glass, amber narrow mouth flat bottom	1.5 ml screwcap vial 32*11.6 mm glass, clear narrow mouth flat bottom write on patch	1.5 ml screwcap vial 32*11.6 mm glass, amber narrow mouth flat bottom write on patch
Packaging	100 pcs / box	100 pcs / box	100 pcs / box	100 pcs / box

			
Article number	120018	120019	120010
Description	0.2 ml screwcap vial 32*11.6 mm glass, clear flat bottom write on patch integrated insert	0.2 ml screwcap vial 32*11.6 mm glass, amber flat bottom write on patch integrated insert	0.3 ml screwcap vial 32*11.6 mm PP, natural flat bottom conical
Packaging			