

# Vial 9 mm screwcap



Article	Description	Packaging
120018	Vial, screwcap, 9 mm, 0.2 ml, 32*11.6 mm, glass, clear, flat bottom, writing patch, integrated insert	100
120019	Vial, screwcap, 9 mm, 0.2 ml, 32*11.6 mm, glass, amber, flat bottom, writing patch, integrated insert	100
120010	Vial, screwcap, 9 mm, 0.3 ml, 32*11.6 mm, PP, slightly transparent, conical, flat bottom	100
120009	Vial, screwcap, 9 mm, 1.5 ml, 32*11.6 mm, glass, clear, flat bottom, NM	100
120023	Vial, screwcap, 9 mm, 1.5 ml, 32*11.6 mm, glass, clear, flat bottom, WM	100
120016	Vial, screwcap, 9 mm, 1.5 ml, 32*11.6 mm, glass, amber, flat bottom	100
120014	Vial, screwcap, 9 mm, 1.5 ml, 32*11.6 mm, glass, clear, flat bottom, writing patch	100
120015	Vial, screwcap, 9 mm, 1.5 ml, 32*11.6 mm, glass, amber, flat bottom, writing patch	100
120024	Vial, screwcap, 9 mm, 1.5 ml, 32*11.6 mm, glass, amber, flat bottom, writing patch, silanized	100
120100	Vial, screwcap, 9 mm, 1.5 ml, 32*11.6 mm, PP, slightly transparent, flat bottom	100



151144



151145



151148



151149



151171

Article	Description	Packaging
151144	Cap, screw, for vial, 2 ml, 9 mm, blue, polypropylene, Sil/PTFE, white/red, packed/100, with hole	1000
151145	Cap, screw, for vial, 2 ml, 9 mm, transp., polypropylene, Sil/PTFE, white/red, packed/100, with hole	1000
151148	Cap, screw, for vial, 2 ml, 9 mm, red, polypropylene, Sil/PTFE, white/red, packed/100, with hole	1000
151149	Cap, screw, for vial, 2 ml, 9 mm, zwart, polypropylene, Sil/PTFE, white/red, packed/100, with hole	1000
151171	Cap, screw, for vial, 2 ml, 9 mm, PP, blue, Sil/PTFE, white/red, 1,0 mm, pre-slit	1000