

Syringe Filters 70 mm

Designed to facilitate the sample preparation of water samples



SPECIFICATIONS

Membrane:	PES (Polyethersulfone)
Pore Size:	0.45 μm
Wettability:	Hydrophilic
Prefilter:	Multilayer Polypropylene
Sterile:	Non sterile
Housing:	Polypropylene
Diameter:	70 mm
Filtration Area:	25 cm ²
Connector in-/ outlet:	8~12 mm stepped hose barb

SPECIAL FEATURES

- High dirt loading capacity
- Low protein binding
- Wide chemical compatibility
- Biologically inert
- Inlet is marked on the device



ISO 9001 ISO 14001

BUREAU VERITAS
Certification



Product number	Description	Packaging
471471	Filter, syringe, PES, 0.45 μm , 70 mm, blue	50 pieces/package
471472	Filter, syringe, PES + prefilter, 0.45 μm , 70 mm, yellow	50 pieces/package